



## Product Data

<b>Product:</b>	Chinese Hamster Hybloc™ DNA
<b>Catalog #:</b>	CHHB
<b>Concentration:</b>	1 mg/mL (A <sub>260</sub> )
<b>Source:</b>	Male Chinese Hamster Tissues
<b>Storage:</b>	-20°C

### **Description:**

Chinese Hamster Hybloc™ DNA is the repetitive sequence fraction of genomic DNA obtained from male Chinese hamster tissues by extracting, shearing, denaturing, and reannealing DNA under conditions that enrich for repetitive sequences. Chinese Hamster Hybloc™ DNA may be effectively substituted for any other competitor DNA (e.g., COT-1 DNA®, salmon sperm DNA, genomic DNA, etc.) to suppress cross-hybridization to Chinese hamster repetitive sequences for *in situ* hybridization, Southern blots, microarray technology, or any other protocol requiring competition to block repetitive DNA from binding to the target. Hybridization with Hybloc™ DNA shows a marked increase in signal intensity and a significant decrease in the amount of background noise associated with repetitive DNA.

### **How Supplied:**

Chinese Hamster Hybloc™ DNAs are supplied in various quantities at 1 mg/mL in 10 mM Tris-HCl (pH 7.4), 1 mM EDTA.

### **Quality Control Assay:**

Determination of concentration from the A<sub>260</sub> and determination of the A<sub>260</sub>/A<sub>280</sub> ratio (1.8 to 2.0); determination of molecular size: ≥80% of Chinese Hamster Hybloc™ DNA will run as a broad band between the 50 bp and 600 bp size range.

### **References:**

- Landegent JE, Jansen in de Wal N, Dirks RW, Bao F, van der Ploeg M. Use of whole cosmid cloned genomic sequences for chromosomal localization by non-radioactive *in situ* hybridization. *Human Genetics* 1987 Dec;77(4):366-70.
- Lengauer C, Riethman H, Cremer T. Painting of human chromosomes with probes generated from hybrid cell lines by PCR with Alu and L1 primers. *Human Genetics* 1990 Nov;86(1):1-6.
- Lichter P, Cremer T, Borden J, Manuelidis L, Ward DC. Delineation of individual human chromosomes in metaphase and interphase cells by *in situ* suppression hybridization using recombinant DNA libraries. *Human Genetics* 1988 Nov;80(3):224-34.
- Lichter P, Tang CJ, Call K, Hermanson G, Evans GA, Housman D, Ward DC. High resolution mapping of human chromosome 11 by *in situ* hybridization with cosmid clones. *Science* 1990 Jan5;247(4938):64-9.
- Scalzi JM, & Hozier JC. Comparative genome mapping: mouse and rat homologies revealed by fluorescence *in situ* hybridization. *Genomics* 1998 Jan1;47(1):44-51.

®COT-1 DNA is a registered trademark of Invitrogen Corporation.